

Install and Setup of Sentinel Medic Patch Driver

After you have installed the DSI software please follow these steps to configure your parallel security key:

- 1) On the Windows desktop, click START button. Go to 'Programs', then 'Rainbow Technologies' group, then the 'Sentinel Medic' group and click on **Sentinel Medic application**.
- 2) A window will pop up. It will ask: "Would you like to install parallel port driver?" Make sure the parallel port driver is selected and click yes/ok.
- 3) Now you should get the message: "Driver installed successfully. No reboot required." If you do not get this message, repeat step 2.
- 4) Now the Sentinel Medic is open. Go to 'Select Key' drop down menu. Click on 'Scribe'. Then click on the large 'Find Keys' button.
- 5) At the bottom you will see the following:
Error and Description sections. Under each of these headings is a box.
*In the Error box, make sure it reads: **0**
*In the Description box make sure it reads:
Sentinel Scribe detected on LPT1
- 6) Close the Sentinel Medic. You will not need to do this again. Run DSI software.



If you don't have a Program group 'Rainbow Technologies' in the Start Menu of Windows, please go to our web site at: <http://www.displaysystemsintl.com/download/utility.htm> and download the Sentinelmedic.exe file. Once the download is complete, run this file. Then proceed to step 1.